

# NGS – variant analysis

Sequencing and alignment





Image from: G. Silva (2016)

# Illumina sequencing


- Sequencing-by-synthesis: 2nd generation sequencing
- Massive throughput: up to  $500 \times 10^9$  bases/run
- Most used platform today

# Illumina sequencing

- 50 – 300 bp
- Paired-end (or single-end)



# Illumina library prep



# fastq


fasta + basequality (fasta + q = fastq)

$$BASEQ = -10\log_{10} \Pr\{base\ is\ wrong\}$$

$$-10\log_{10} (0.01) = 20$$

$$-10\log_{10} (0.1) = 10$$

$$-10\log_{10} (0.5) = 3$$



# fastq

reads.fastq

```
@D00283R:66:CC611ANXX:4:2311:2596:2330 1:N:0:TCCGGAG
ACTCTACGCTCAATAAAGATTTCTGATACGGCTCCTGAAATGCAGAATGAGT
+
B/<<<B<FFFFFFFFFBBFFFBFFFBFFFF/FFFFFFFF/BFFFBFFF
```


title, starts with @

nucleotide sequence

optional description

base quality


# Base quality (phred)





S - Sanger Phred+33, raw reads typically (0, 40)  
X - Solexa Solexa+64, raw reads typically (-5, 40)  
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)  
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)  
with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)  
(Note: See discussion above).  
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)




# Quiz Question 6






in phase



out of phase


# Illumina - limitations

- Bridge amplification
- Lengths are limited by out-of-phase of signal



# Long reads (3rd generation)

- Crux: maximizing signal from a single-molecule base read-out
- Single molecule, so no out-of-phase signal
- Two frequently used platforms:
  - PacBio SMRT sequencing
  - Oxford Nanopore Technology



# Long reads


- More error -> difficulties for variant analysis
- But:
  - PacBio CCS: high baseQ + no bias
  - Long reads can have higher mapping qualities
  - Long reads improve haplotyping


# What to sequence?

- Whole genome/metagenome
- Reduced representation:
  - Bait capture/whole exome sequencing
  - Restriction Enzyme based (e.g. RAD seq)
  - Amplicon sequencing
  - RNA-seq



# Read alignment (phred)






# Software

- Basic alignment:
  - bowtie2
  - bwa-mem
- Long reads:
  - minimap2


# Mapping quality



$$MAPQ = -10\log_{10} \Pr\{\text{mapping position is wrong}\}$$

$$-10\log_{10} (0.01) = 20$$

$$-10\log_{10} (0.5) = 3$$



# Quiz Question 7

# sam

## sequence alignment format



# sam header


```
@HD      VN:1.0  SO:coordinate
@SQ      SN:U00096.3      LN:4641652
@PG      ID:bowtie2      PN:bowtie2      VN:2.4.1      CL: bowtie2-
align-s --wrapper basic-0 -x ref.fasta -1 reads_1.fastq -2
reads_2.fastq"
```

<b>SAM column</b>	<b>example</b>
read name	SRR519926.5
flag	89
reference	chr20
start position	61
mapping quality	42
CIGAR string	150M
reference name mate is mapped	=
start position mate	476
fragment length	515
sequence	CATCACCATTCCCAC
base quality	@>4:4C@89+&9CC@
optional	AS:i:-2
optional	XN:i:0



# Quiz Question 8

# Marking duplicates



# Marking duplicates

- Variant calling: each read is an independent observation of the genome
- Duplicates (can) have the same molecular origin -> not independent
- Removing duplicates probably doesn't have a big effect on variant analysis

Ebbert MTW et al. (2016) Evaluating the necessity of PCR duplicate removal from next-generation sequencing data and a comparison of approaches. BMC Bioinformatics.

# Unique Molecular Identifiers

- UMI added before PCR reaction
- Detect PCR duplicates and PCR errors