

NGS – variant analysis

Sequencing and alignment




Image from: G. Silva (2016)

Illumina sequencing


- Sequencing-by-synthesis: 2nd generation sequencing
- Massive throughput: up to 500×10^9 bases/run
- Most used platform today

illumina[®]



Illumina sequencing

- 50 – 300 bp
- Paired-end (or single-end)



Illumina library prep

shear DNA



Ligate adapters



Barcode + p5/p7 sites



PCR: 8-16 cycles



Sequencing



fastq

fasta + basequality (fasta + q = fastq)

$$BASEQ = -10\log_{10} \Pr\{base\ is\ wrong\}$$

$$-10\log_{10} (0.01) = 20$$

$$-10\log_{10} (0.1) = 10$$

$$-10\log_{10} (0.5) = 3$$

fastq

reads.fastq

```
@D00283R:66:CC611ANXX:4:2311:2596:2330 1:N:0:TCCGGAG
ACTCTACGCTCAATAAAGATTTCTGATACGGCTCCTGAAATGCAGAATGAGT
+
B/<<<B<FFFFFFFFFBBFFFBFFFBFFFF/FFFFFFFF/BFFFBFFF
```


title, starts with @

nucleotide sequence

optional description

base quality

Base quality (phred)





S - Sanger Phred+33, raw reads typically (0, 40)
X - Solexa Solexa+64, raw reads typically (-5, 40)
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)
with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
(Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)

Quiz Question 4

With `grep -c "^>"` it is easy to count the number of sequences in a fasta file. Why doesn't `grep -c "^@"` work for fastq sequences?



- A. @ is a special character for regular expressions
- B. The @ can also occur elsewhere in a fastq file
- C. There is no specific title start symbol for fastq files





Illumina - limitations


- Bridge amplification
- Lengths are limited by out-of-phase of signal



Long reads (3rd generation)

- Crux: maximizing signal from a single-molecule base read-out
- Single molecule, so no out-of-phase signal
- Two frequently used platforms:
 - PacBio SMRT sequencing
 - Oxford Nanopore Technology





Question 4

Why are INDELs more difficult to detect from alignments compared to SNPs?


- A. Because there is no base quality
- B. Because it is difficult to know the correct local alignment


Long reads

- More error -> difficulties for variant analysis
- But:
 - PacBio CCS: high baseQ + no bias
 - Long reads can have higher mapping qualities
 - Long reads improve haplotyping



Read alignment (phred)






Software

- Basic alignment:
 - bowtie2
 - bwa-mem
- Long reads:
 - minimap2

Mapping quality



$$MAPQ = -10\log_{10} \Pr\{\text{mapping position is wrong}\}$$

$$-10\log_{10} (0.01) = 20$$

$$-10\log_{10} (0.5) = 3$$

sam

sequence alignment format



sam header

```
@HD      VN:1.0  SO:coordinate
@SQ      SN:U00096.3      LN:4641652
@PG      ID:bowtie2      PN:bowtie2      VN:2.4.1
CL: "/opt/miniconda3/envs/ngs/bin/bowtie2-align-s \
--wrapper basic-0 \
-x /home/ubuntu/ecoli/ref_genome//ecoli-strK12-MG1655.fasta \
-1 /home/ubuntu/ecoli/trimmed_data/paired_trimmed_SRR519926_1.fastq \
-2 / home/ubuntu/ecoli/trimmed_data/paired_trimmed_SRR519926_2.fastq"
```


SAM column	example
read name	SRR519926.5
flag	89
reference	U00096.3
start position	61
mapping quality	42
CIGAR string	214M
reference name mate is mapped	=
start position mate	476
fragment length	515
sequence	CATCACCATTCCCAC
base quality	@>4:4C@89+&9CC@
optional	AS:i:-2
optional	XN:i:0

Question 5

Can you technically regenerate the fastq file out of the SAM file? And can you regenerate the reference sequence (fasta) file from the SAM file?

- A. Only the fastq file
- B. Only the fasta file
- C. Both files
- D. None of those

Marking duplicates



Marking duplicates

- Variant calling: each read is an independent observation of the genome
- Duplicates (can) have the same molecular origin -> not independent
- Removing duplicates probably doesn't have a big effect on variant analysis

Ebbert MTW et al. (2016) Evaluating the necessity of PCR duplicate removal from next-generation sequencing data and a comparison of approaches. BMC Bioinformatics.

Unique Molecular Identifiers

- UMI added before PCR reaction
- Detect PCR duplicates and PCR errors